



HEALTH HOLDING

HAFER ALBATIN HEALTH
CLUSTER
MATERNITY AND
CHILDREN HOSPITAL

Department:	Anesthesia Care		
Document:	Departmental Policy and Procedure		
Title:	Cleaning/Disinfecting of Anesthesia Machine		
Applies To:	All Anesthesia Staff, Anesthesia Technician and Nursing Staff in Operating Room		
Preparation Date:	January 05, 2025	Index No:	AN-DPP-017
Approval Date:	January 19, 2025	Version :	2
Effective Date:	February 19, 2025	Replacement No.:	AN-DPP-017(1)
Review Date:	February 19, 2028	No. of Pages:	03

1. PURPOSE:

- 1.1 To establish quality standard in disinfecting and cleaning procedure for anesthesia machine to prevent infection transmission to patient.

2. DEFINITIONS:

- 2.1 This is the cleaning and disinfection process and manual for anesthesia machine.

3. POLICY:

- 3.1 **Anesthesia Technicians** – are responsible for the cleaning and disinfecting anesthesia machines, anesthesia breathing circuits and surfaces and monitoring according to the next procedure.

4. PROCEDURE:

4.1 The Anesthesia Machine

Although there is no direct contact between anesthesia machine controls and the patient, microorganisms can be transferred between the machine and patient by the healthcare provider.

- 4.1.1 Cleaning and disinfecting the anesthesia machine.
- 4.1.2 Pasteurizing or autoclaving of valves.
- 4.1.3 Disassembling and disinfecting adjustable pressure-limiting valves.

4.2 Anesthesia Machine Surfaces and Carts

Clean, then spray or wipe anesthesia machine surfaces and knobs with an appropriate germicide between cases and at the end of each day.

- 4.2.1 Take protective measures to prevent materials stored on the anesthesia machine from becoming inadvertently contaminated by airborne debris (e.g., blood).
- 4.2.2 Remove equipment from drawers, clean and disinfect drawers regularly.
- 4.2.3 Place a clean covering on the top of the anesthesia cart at the beginning of each case.
- 4.2.4 Wipe small surfaces with 70 percent isopropyl alcohol to reduce bacterial contamination.
- 4.2.5 Clean carbon dioxide and soda lime absorbers when the absorber is changed and remove debris from the screens.

4.3 Anesthesia Breathing System

Filters Breathing System – filters are single – use items that are assessed according to their bacterial filtration efficiency (BFE) and viral filtration efficiency (VFE). The efficacy of filtration for bacterial contaminants is higher than for viral particles. Aside from patients with an active Mycobacterium Tuberculosis Infection, no recommendation is made for the routine use of breathing system filters due to inconclusive data demonstrating their efficacy in reducing the risk of patient infection. However, when a patient with a respiratory infection must be given inhalational anesthesia, a filter should be used.

- 4.3.1 Practitioners may choose to place a high-efficiency filter on the inspiratory limb of the breathing circuit to protect the patient from the anesthesia machine, and to place a high-efficiency filter in

the expiratory limb to protect the anesthesia machine from the patient. • Filters may be interposed between the endotracheal tube and the Y-piece.

- 4.3.2 Use circuit filters and follow-up with post-anesthesia machine disinfection after caring for patients with known pulmonary infection or trauma.

4.4 **Carbon Dioxide Absorbers**

- 4.4.1 Cleaning, and sterilization of carbon dioxide absorbers :

- 4.4.1.1 Clean canisters when the absorbent is changed and carefully remove debris from the screens.
- 4.4.1.2 Discard disposable plastic canisters.
- 4.4.1.3 Bellows, unidirectional valves, and carbon dioxide absorbers should be cleaned and disinfected periodically.

4.5 **Circuits Anesthesia**

- 4.5.1 Circuits may be manufactured as either single patient use items or multiple patient use items (provided that a new breathing system filter is placed between the Y-piece and endotracheal tube after sterilization or high level disinfection).
- 4.5.2 Anesthesia professionals should pay close attention to anesthesia circuit product labelling.
 - 4.5.2.1 At a minimum, provide high-level disinfection for multiple-patient use breathing circuits.
 - 4.5.2.2 The outer surface of the circuit can become easily contaminated when the system is not changed between patients and therefore should be disinfected between each use.
 - 4.5.2.3 End- tidal carbon dioxide tubing should be changed between patients.
 - 4.5.2.4 Following anesthesia care of a patient with pulmonary infection or trauma, disinfection of the internal and respiratory system anesthesia machine components is mandatory.
- 4.5.3 Heat and Moisture Exchangers: Heat and moisture exchangers alone are not effective in decreasing the transmission of microorganisms to the anesthesia breathing system.

5. **MATERIALS AND EQUIPMENT:**

- 5.1 The approved cleaning materials from IPC department.

6. **RESPONSIBILITIES:**

- 6.1 Nurses
- 6.2 Anesthesia Technician









7. **APPENDICES:**

N/A

8. **REFERENCES:**

- 8.1 American Society of Anesthesiology.

9. APPROVALS:

	Name	Title	Signature	Date
Prepared by:	Ms. Naimah Naif Al Salem	Head Nurse of the Operating Room		January 05, 2025
Prepared by:	Mr. Hamed Matar Alanazi	Head of Anesthesia Technician		January 05, 2025
Prepared by:	Dr. Abdelghani Ibrahim	Head of the Anesthesia Department		January 05, 2025
Reviewed by:	Mr. Sabah Turayhib Al Harbi	Director of Nursing		January 08, 2025
Reviewed by:	Ms. Awatif Hamoud Al Harbi	IPCD Director		January 09, 2025
Reviewed by:	Mr. Abdulelah Ayed Al Mutairi	QM&PS Director		January 09, 2025
Reviewed by:	Dr. Tamer Mohamed Naguib	Medical Director		January 12, 2025
Approved by:	Mr. Fahad Hazam Al - Shammari	Hospital Director		January 19, 2025